

Take a closer look at the Dräger Fabius Tiro

In the modern hospital environment, time is always precious. But there are also many places in a hospital where space comes at a premium. Places like induction rooms, emergency departments or smaller operating rooms. Spaces like these present particular challenges- such as how to maintain the best possible patient access and still provide high-quality care. It was for environments like these that we designed the Dräger Fabius Tiro

The Fabius Tiro combines everything you need to provide high-quality ventilation in small spaces. Its compact design makes it a natural choice for applications where you need to get the most of the available room without having to compromise on features or flexibility. With a wide array of options to choose from, the Fabius Tiro can be customized to suit your needs and even be upgraded when those needs change.



Performance Features:

- Compact, space-saving design
- Highly maneuverable cart
- Optional wall mount
- Standardized Dräger user interface for easy and intuitive operation
- All major ventilation modes
- High-contrast display
- Intelligent safety features for enhanced patient protection
- CLIC absorber compatible

» HIGH-PERFORMANCE VENTILATION

- Powerful and accurate E-vent piston ventilation which requires no drive gas and delivers ICU-like ventilation performance.
- Fresh gas decoupling which maintains constant tidal volumes regardless of fresh gas flow changes.
- Compliance compensation and backup ventilations for PS modes let you ventilate with enhanced safety, confidence, and control.