

# SKYTRON AURORA SERIES

*Fully Refurbished*



**Heartland**  
Medical Sales & Services, LLC

HeartlandMedical.com

The Skytron Aurora LED surgical lights provide an unmatched color-correct illumination with advanced bright and cool light and a selectable color temperature. Based on preference, surgeons have a choice of 4,000K (soft white) or 4,500K (bright white) for every procedure. Focusable LEDs extend illumination range while retaining maximum intensity at any distance. Common applications include: General, Orthopedic, Plastics, Vascular, Cardiovascular, OB/GYN, Trauma, Neurosurgery, and Urology.



## FEATURES:

- Direct surgeon control of intensity & focus
- Brighter intensity, deep-view illumination
- Cool beam temperature with less heat
- Shadow-free optics
- Ultra energy-efficient & low maintenance
- Eco-friendly LEDs last over a decade
- Precise, 360° positioning at all points of articulation, without drift
- Advanced architecture
- HD camera ready
- Color correct temperature of 4,000K or 4,500K can be determined case-by-case or during the surgical case based upon the surgeon's needs and preference
- The AUR 5 has 5 LED pods within its head to provide a light up to 135,000 lux
- The Aurora II uses 3rd generation LED technology to provide optimally bright, cool and amalgamated white light
- Dynamically focusable satellite LEDs combined with high-intensity fixed-focus Centra LEDs deliver superior, deep cavity illumination, without shadows
- Advanced Vertically Segmented Reflector Design (VSRD) optics greatly enhance the lighting performance of each LED



HEARTLAND MEDICAL SALES & SERVICES, LLC

2601 HOLLOWAY ROAD LOUISVILLE, KY 40299 | 1-502-671-1014 | 877-329-7775

[www.HeartlandMedical.com](http://www.HeartlandMedical.com)